



# Trillium 2000™ Holter Recorder

## Add Holter to Your Practice

Now you can quickly and easily add on-site Holter monitoring capabilities to your practice with our Trillium 2000 Holter Recorders. They are extremely rugged and come with a two year warranty and unlimited toll free technical support from our Holter knowledgeable staff.

Trillium 2000 Holter Recorders can be combined with our Trillium Silver™, Trillium Gold™ or Trillium Platinum™ Holter Analysis Software to form a complete Holter Monitoring System that is easy to use. All of our compatible software provides fast and accurate Holter analysis and reporting. Please see our Trillium Silver, Trillium Gold, and Trillium Platinum data sheets for more details

## Qualify for Maximum Reimbursement

Trillium 2000 based Holter Monitoring Systems can be used to provide Holter diagnostic procedures that qualify for maximum reimbursement (CPT Code 93230).

## Enhance Patient Retention

The diagnostic quality reports produced on-site by our Trillium 2000 based Holter Monitoring Systems can help you to determine if a referral to a specialist is actually required. Reducing the number of unnecessary referrals enhances patient retention while improving patient service through faster turnaround times, fewer visits, and broader point of care offerings. Improving patient retention rates can help you to maintain and build your practice while increasing your bottom line.

## Discover the Value Leader in Holter

Our Holter monitoring products are surprisingly affordable and provide the best long term value in Holter monitoring technology. Our systems are easily and economically expandable and are fully upgradable. By choosing Forest Medical, you can be sure that you've made the right investment in Holter monitoring technology now and for the future.



## Trillium 2000 Holter Recorder Specifications

- Channels: 2 or 3
- Recording Time: 24 Hours
- Recording: Continuous (Full Disclosure)
- Memory: Non-volatile Flash-Card
- Resolution: 8 – 12 bit sampling
- Sample Rate: 256 samples per second per channel
- Frequency Response: 0.05 Hz – 100 Hz
- Gain Settings: Automatic per channel
- Power Requirements: (3) AA Alkaline IEC-LR6
- Overall size: 3.5 x 6.4 x 0.9 inches (88.9 x 162.6 x 22.8 mm)
- Weight Including Batteries: 10.75 ounces (305 grams)
- Operating Temperature: 32° F to 113° F (0° C to 45° C)
- Non-operating Temperature: -4° F to 149° F (-20° C to 65° C)
- Operating Humidity: 10% to 95% (non-condensing)
- Non-operating Humidity: 5% to 95% (non-condensing)
- Patient Cable: 5 or 7 Lead
- Signal Verification: Multi-color LEDs and Test Strips
- Three button operation with 16 character LCD display
- Time and Parameter Settings: 5 Year Lithium battery backup
- Meets the ANSI/AAMI standard for ambulatory electrocardiographs (EC38)
- Unlimited Toll Free Technical Support

Notice: In the United States, Federal law restricts this device to sale by or on the order of a physician.

## Call Forest Medical

Call us today to get started with Holter monitoring at your practice.

Widely Accepted  
Diagnostic Procedure

Rugged and Reliable

24 Hour Continuous  
Recording

Choice of PC Based  
Software for Analysis

Qualifies for Maximum  
Reimbursement

Fast Return  
on Investment

Thousands in Use

Easy to Buy via  
COD, Credit Card, or Lease

Quick Delivery

2 Year Warranty

Turnkey Systems Available

Forest Medical, LLC • 6700 Old Collamer Road • East Syracuse, NY 13057  
Phone: 800-844-2037 • 315-434-9000 • Fax: 315-432-8064  
sales@forestmedical.com • www.forestmedical.com