



Trillium 5000™ Holter Recorder

Advanced Technology

Compact Recorder

Pacemaker Detection

24 or 48 Hour
Continuous Recording

Choice of PC Based
Software for Analysis

Fast Return
on Investment

Customizable
Leasing Programs

1 Year Warranty

Turnkey Systems Available

Choose Compact Versatility

Our compact and versatile Trillium 5000 Holter Recorders provide 2 or 3 channel, 24 or 48 hour, recording capability. They provide complete support for pediatric through adult applications including pacemaker detection. Trillium 5000 Holter Recorders provide superior performance and come with a one year warranty and unlimited toll free technical support from our Holter knowledgeable staff.

Trillium 5000 Holter Recorders can be combined with our Trillium Silver™, Trillium Gold™, or Trillium Platinum™ Holter Analysis Software to form a complete and comprehensive Holter Monitoring System that is both powerful and easy to use.

Provide the Best Possible Care

Our Trillium Silver, Trillium Gold, and Trillium Platinum Holter Analysis Software products provide flexible and accurate analysis, expanded functionality, powerful editing capabilities, and configurable and editable reports to give you the information you need to provide the best possible patient care. And, with a Trillium 5000 based Holter Monitoring System you can choose the level of software that best meets the requirements of your practice. Please see our Trillium Silver, Trillium Gold, and Trillium Platinum Holter Analysis Software data sheets for more details.

Discover the Value Leader in Holter

Forest Medical Holter Monitoring Systems are surprisingly affordable and provide the best long term value in Holter monitoring technology. Our Holter Monitoring Systems are easily and economically expandable with multiple recorders and workstations and are fully upgradable as new developments arise. By choosing Forest Medical, you can be sure that you've made the right investment in Holter monitoring technology now and for the future.

Call Forest Medical

Call us today for a no obligation analysis of your Holter monitoring requirements.

Forest Medical, LLC • 6700 Old Collamer Road • East Syracuse, NY 13057
Phone: 800-844-2037 • 315-434-9000 • Fax: 315-432-8064
sales@forestmedical.com • www.forestmedical.com



Trillium 5000 Holter Recorder Specifications

- Channels: 2 or 3
- Recording Time: 24 or 48 Hours
- Recording: Continuous (Full Disclosure)
- Memory: Non-volatile Flash-Card
- Resolution: 8 bit sampling
- Sample Rate: 256 samples per second per channel
- Frequency Response: 0.05 Hz – 60 Hz
- Gain Settings: 1/2X, 1X, and 2X
- Power Requirements: (1) AA Alkaline IEC-LR6
- Overall size: 2.78 x 3.75 x 0.78 inches (88.9 x 162.6 x 24.1 mm)
- Weight Including Battery: 4 ounces (112 grams)
- Operating Temperature: 32° F to 113° F (0° C to 45° C)
- Non-operating Temperature: -4° F to 149° F (-20° C to 65° C)
- Operating Humidity: 10% to 95% (non-condensing)
- Non-operating Humidity: 5% to 95% (non-condensing)
- Patient Cable: 5 or 7 Lead
- Signal Verification: LCD Display - ECG Display on LCD
- Five button keypad operation with LCD Display
- Meets the ANSI/AAMI standard for ambulatory electrocardiographs (EC38) including pacemaker and use on infants
- Unlimited Toll Free Technical Support

Notice: In the United States, Federal law restricts this device to sale by or on the order of a physician.